

Eleiko Prestera Hip Thrust Pad

Galvanized/PUR

Designed for the Prestera System, the Hip Thrust Pad attaches between two uprights creating a safe, dedicated space for lifters to perform hip thrust and glute exercises. The attachment is easily placed between the uprights and securely fastened with Eleiko Magnetic Pins.



Funktioner

Nem at bruge og opbevare

Hold puden i den ønskede højde, og fastgør den med magnetiske stifter. Hvis puden opbevares adskilt fra raket, kan stifterne opbevares i tilbehøret. Ellers kan de monteres mellem stolperne, hvor de ikke blokerer for noget.

Træningsvariation

Velegnet til hoftepresøvelser, der udføres med vægte, bånd eller anden belastning. Puden kan også bruges til bulgarske split-squats og er robust nok til step-ups.

Kvalitetskonstruktion

Stærke stålrør og rustfrit galvaniseret stål giver langvarig ydeevne udendørs, og PUR-puden af polyurethan er nem at rengøre og vedligeholde.

Kompatibilitet

Designet specielt til Prestera-serien.

Specifikationer

Artikelnummer

BD-1-368

Længde

1231,2 mm / 48.47 in.

Bredde

528 mm / 20.79 in.

Højde

150 mm / 5.91 in.

Vægt

10,25 kg / 22.6 lbs

Garanti

1-10 år

