

# Eleiko Prestera Hip Thrust Pad

Black/PUR

Conçu pour s'intégrer à la gamme Prestera, le module pour Hip Thrust se fixe entre deux montants pour créer un espace sécurisé et réservé aux hip thrust et autres mouvements musclant les fessiers. Il est facile à positionner et à fixer entre deux poteaux, à l'aide des goupilles magnétiques Eleiko.



## Fonctionnalités

### Pratique

Simply hold the pad at desired heights and secure it with Eleiko Magnetic Pins. If storing the unit separately from the rack, the pins can be stored in the attachment. Alternatively, mount it out of the way between uprights.

### Fabrication de Qualité

Sturdy steel tubing with durable powder coat finish and polyurethane PUR pad is easy to clean and maintain.

### Exercices variés

Appropriate for hip thrust exercises performed with barbells, bands, or other loads. It can also be used for elevated rear foot split squats and is sturdy enough for step-ups.

### Compatibilité

Spécialement conçu pour la gamme Prestera.

## Caractéristiques

### Code Article

BD-1-367

### Longueur

1291 mm / 50.83 in.

### Largeur

400 mm / 15.75 in.

### Hauteur

168 mm / 6.61 in.

### Poids

13,7 kg / 30.2 lbs

### Garantie

6 mois à 10 ans

